

PHILIPS

Supplies catalog 2025





Directory



Anesthesia gas

Anesthesia gas monitoring

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
451261001391	M1662A	View >	Anesthetic calibration gas (use with M1659A) Volume of gas: 6 liters Shelf life: 2 years from the date of manufacture	1 cylinder	Anesthesia gas module: G7m 866173	•	•
989803100361	13907A	View >	Calibration tube assembly	1 assembly	Anesthesia gas analyzer: G7m 866173, Microstream™ ext.: 862393, 865377	•	•
989803100841	15210-64010	View >	Gas: Cal 1 cylinders for tcpCO ₂	6 cylinders	Microstream™ ext.: 862393, 865377	•	•
989803100851	15210-64020	View >	Gas: Cal 2 cylinders for tcpCO ₂	6 cylinders	Microstream™ ext.: 862393, 865377	•	•
989803104681	M1659A	View >	Calibration bag (use with M1662A)	1 bag	Anesthesia gas module: G7m 866173	•	•
989803191081		View >	Water trap	12 water traps	G7m: 866173	•	•



451261001391
M1662A



989803100361
13907A



989803100841
15210-64010



989803100851
15210-64020



989803104681
M1659A



989803191081



Spirometry CO₂

IntelliVue Spirometry CO₂ sensors

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803144661	M2785A	Flow sensor, adult/pediatric, single-patient use ET: > 4.0 mm (0.16 in)	10 sensors	Spirometry module: 865021	•	•
989803144671	M2786A	Flow sensor, infant/neonatal, single-patient use ET: < 4.0 mm (0.16 in)	10 sensors	Spirometry module: 865021	•	•



989803144661
M2785A



989803144671
M2786A

Sidestream CO₂

LoFlo Sidestream CO₂ sensor

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803144591	M2741A	Sidestream CO ₂ sensor	1 sensor	Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; XPER Flex 860335	•	•



989803144591
M2741A


Sidestream CO₂

LoFlo Sidestream CO₂ intubated

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206711	n/a	Airway Adapter - LoFlo - ET >4.0 mm	10 sets	Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; XPER Flex 860335	•	•
989803206721	n/a	Airway Adapter - LoFlo - ET <4.0 mm	10 sets		•	•



989803206711



989803206721

Sidestream CO₂

LoFlo Sidestream CO₂ intubated with humidification

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206731	n/a	Airway Adapter Set H - ET >4.0 mm	10 sets	Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; XPER Flex 860335	•	•
989803206741	n/a	Airway Adapter Set H - ET <4.0 mm	10 sets		•	•



989803206731

989803206741

Sidestream CO₂

LoFlo Sidestream CO₂ non-intubated with O₂ delivery

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206641	n/a	Cannula - LoFlo - CO ₂ /O ₂ Nasal - Adult	10 cannulas	Sidestream LoFlo sensor: M2741A, 865040; MP5: 865024	•	•
989803206651	n/a	Cannula - LoFlo - CO ₂ /O ₂ Nasal - Pediatric	10 cannulas		•	•
989803206661	n/a	Cannula - LoFlo - CO ₂ /O ₂ Nasal - Infant	10 cannulas		•	•
989803206691	n/a	Cannula - LoFlo - CO ₂ /O ₂ Oral/Nasal-Adult	10 cannulas		•	•
989803206701	n/a	Cannula - LoFlo - CO ₂ /O ₂ Oral/Nasal-Pediatric	10 cannulas		•	•



989803206641




989803206651



989803206661



989803206691



989803206701

Sidestream CO₂

LoFlo Sidestream CO₂ non-intubated without O₂ delivery

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206611	n/a	Cannula - LoFlo - CO ₂ Nasal - Adult	10 cannulas	Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187, X2: 865039; MP2: 865040; MP5: 865024	•	•
989803206621	n/a	Cannula - LoFlo - CO ₂ Nasal - Pediatric	10 cannulas		•	•
989803206631	n/a	Cannula - LoFlo - CO ₂ Nasal - Infant	10 cannulas		•	•
989803206671	n/a	Cannula - LoFlo - CO ₂ Oral/Nasal - Adult	10 cannulas		•	•
989803206681	n/a	Cannula - LoFlo - CO ₂ Oral/Nasal-Pediatric	10 cannulas		•	•



989803206611



989803206621



989803206631



989803206671



989803206681

Sidestream CO₂ LoFlo Sidestream CO₂ supplies

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206751	n/a	Sample Line-LoFlo - w/Male Luer	10 lines	Sidestream LoFlo sensor: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; XPER Flex 860335	•	•
989803206761	n/a	Sample Line-LoFlo - w/Male Luer-w/Nafion			•	•
989803199061	n/a	LoFlo extension line, length: 1.52 m (60in)	25 lines	Sidestream LoFlo sensor:: M2741A with Capnography ext.: 862187; MP2: 865040; MP5: 865024; X2: 865039; X3: 867030; MX100: 867033; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; XPER Flex 860335		



989803206751



989803206761



989803206631

Sidestream CO₂

CapnoTrak sidestream CO₂ supplies


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803199011	n/a	Water filter assembly	10 pieces	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172	•	
989803199021	n/a	CO ₂ sampling extension line	25 lines		•	
989803199031	n/a	O ₂ delivery extension line	25 lines		•	
989803199041	n/a	Dehumidification tubing	10 tubes		•	
989803199061	n/a	LoFlo extension line, length: 1.52 m (60 in)	25 lines		•	



Sidestream CO₂

CapnoTrak sidestream CO₂ intubated, short-term

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803198991	n/a	CapnoTrak Sidestream CO ₂ airway adapter set, intubated, short-term ET: > 4.0 mm (0.16 in)	25 sets	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172	•	
989803199001	n/a	CapnoTrak Sidestream CO ₂ airway adapter set, intubated, short-term ET: < 4.0 mm (0.16 in)	25 sets		•	



989803198991



989803199001



Sidestream CO₂

CapnoTrak Sidestream CO₂, non-intubated, short-term, with O₂ delivery

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803198921	n/a	CO ₂ / O ₂ nasal cannula, large	25 cannulas	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172	•	
989803198931	n/a	CO ₂ / O ₂ nasal cannula, medium	25 cannulas		•	
989803198941	n/a	CO ₂ / O ₂ nasal cannula, small	25 cannulas		•	
989803198971	n/a	CO ₂ / O ₂ oral-nasal cannula, large	25 cannulas		•	
989803198981	n/a	CO ₂ / O ₂ oral-nasal cannula, medium	25 cannulas		•	



989803198921



989803198931



989803198941



989803198971



989803198981

Sidestream CO₂

CapnoTrak Sidestream CO₂, non-intubated, short-term, without O₂ delivery

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803198891	n/a	CO ₂ / O ₂ nasal cannula, large	25 cannulas	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; Heartstart Intrepid: 867172	•	
989803198901	n/a	CO ₂ / O ₂ nasal cannula, medium	25 cannulas		•	
989803198911	n/a	CO ₂ / O ₂ nasal cannula, small	25 cannulas		•	
989803198951	n/a	CO ₂ / O ₂ oral-nasal cannula, large	25 cannulas		•	
989803198961	n/a	CO ₂ / O ₂ oral-nasal cannula, medium	25 cannulas		•	




989803198891



989803198901



989803198911



989803198951



989803198961

Short term use: procedure/emergency

Microstream™ Advance sampling lines

Ordering code	Image	Description		Packaging	Use with	CE	FDA
Non-intubated patient							
989803204571	View	Oral/nasal sampling line; samples from both the mouth and the nose; integrated O ₂ tubing for supplemental oxygen delivery	Adult; O ₂ male connector; standard: 2 m (6.5 ft)	25 per box	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; CM (Brazil only): 863320; 863321; 863322; Efficia CM Microstream™ Module: 863335	•	•
989803204371	View		Adult; O ₂ male connector; long: 4 m (13 ft)	25 per box		•	•
989803204551	View		Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204361	View		Adult; O ₂ tubing female connector; long: 4 m (13 ft)	25 per box		•	•
989803204561	View	Oral/nasal sampling line; samples from both the mouth and the nose; available with or without O ₂ delivery	Pediatric; no O ₂ delivery; standard: 2 m (6.5 ft) no O ₂ delivery; standard: 2 m (6.5 ft)	25 per box	IntelliVue Microstream™ Extensions 867041, M3015A, M3015B, MP5/MP5SC: M8105A, M8105AS PM6000: 867311, 867313, 867315	•	•
	View		Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204351	View		Pediatric; O ₂ tubing female connector; long: 4 m (13 ft)	25 per box		•	•
989803204481	View	Nasal sampling line; nasal sampling only; integrated O ₂ tubing for supplemental oxygen delivery	Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204491	View		Adult; O ₂ tubing female connector; long: 4 m (13 ft)	25 per box		•	•
989803204501	View		Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204611	View	Nasal sampling line; nasal sampling only; no O ₂ delivery	Adult; standard: 2 m (6.5 ft)	25 per box		•	•
989803204441	View		Infant/neonate; standard: 2 m (6.5 ft)	25 per box		•	•
989803204601	View		Pediatric; standard: 2 m (6.5 ft)	25 per box		•	•
989803204451	View	60fr bite block sampling line; samples from both the mouth and the nose; integrated O ₂ tubing for supplemental oxygen delivery; O ₂ delivery system can administer oxygen flow up to 10 L/min	O ₂ male connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204461	View		O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204471	View		O ₂ tubing female connector; long: 4 m (13 ft)	25 per box		•	•
Intubated patient							
989803204511	View	Intubated sampling line; for use with ETT size > 4.5 mm	Adult/pediatric; standard: 2 m (6.5 ft)	25 per box		•	•
989803204321	View		Adult/pediatric; long: 4 m (13 ft)	25 per box		•	•



989803204571



989803204371
989803204551
989803204361




989803204561
989803204541



989803204481
989803204491
989803204501



989803204611
989803204441
989803204601



989803204451
989803204461
989803204471



989803204511
989803204321

Extended duration Microstream™ Advance sampling lines

Ordering code	Image	Description	Packaging	Use with	CE	FDA	
Non-intubated patient							
989803204381	View	Oral/nasal sampling line; samples from both the mouth and the nose; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery	Adult; O ₂ male connector; standard: 2 m (6.5 ft)	25 per box	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; CM (Brazil only): 863320; 863321; 863322; Efficia CM Microstream™ Module: 863335 IntelliVue Microstream™ Extensions 867041, M3015A, M3015B, MP5/MP5SC: M8105A, M8105AS PM6000: 867311, 867313, 867315	•	•
989803204391	View		Adult; O ₂ male connector; long: 4 m (13 ft)	25 per box		•	•
989803204401	View		Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204411	View		Adult; O ₂ tubing female connector; long: 4 m (13 ft)	25 per box		•	•
989803204621	View	Nasal sampling line; nasal sampling only; drier moisture reduction technology; integrated O ₂ tubing for supplemental oxygen delivery available	Pediatric; no O ₂ delivery; standard: 2 m (6.5 ft) no O ₂ delivery; standard: 2 m (6.5 ft)	25 per box		•	•
989803204581	View		moisture reduction technology; integrated O ₂ tubing	25 per box		•	•
989803204591	View		for supplemental oxygen delivery available	25 per box		•	•
989803204631	View		Adult; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box		•	•
989803204421	View		Adult; O ₂ tubing female connector; long: 4 m (13 ft)	25 per box	•	•	
989803204431	View		Pediatric; O ₂ tubing female connector; standard: 2 m (6.5 ft)	25 per box	•	•	
Intubated patient							
989803204521	View	Intubated sampling line; for use with ETT > 4.5 mm; drier moisture reduction technology		25 per box	•	•	
989803204331	View			25 per box	•	•	
989803204531	View			25 per box	•	•	
989803204341	View			25 per box	•	•	
989803204301	View	Intubated sampling line; for use with ETT ≤ 4.5 mm; longer drier moisture reduction technology for high ambient humidity		25 per box	•	•	
989803204311	View			25 per box	•	•	



989803204321



989803204381
989803204391



989803204401




989803204411




989803204621
989803204581
989803204591




989803204631




989803204521
989803204331



989803204531
989803204341



989803204301



989803204311

Mainstream CO₂ solution

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803142651	M2501A	View >	Mainstream CO, sensor	1 sensor	Capnography ext.: 862187, 867040; MP2: 865040; MP5: 865024, 865322; X2: 865039; X3: 867030; MX100: 867033; MMX 867036; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199; Xper Flex: 860338;	•	•
989803142661	M2533A	View >	Single-patient use adult airway adapter ET: > 4.0 mm (0.16 in)	10 adaptors		•	•
989803142671	M2536A	View >	Single-patient use infant airway adapter ET: <4.0 mm (0.16 in) Deadspace: < 1 cc			•	•
989803206591	n/a	View >	Airway Adapter-Reusable-Adult/Pediatric	1 adaptor		•	•
989803206601	n/a	View >	Airway Adapter - Reusable - Neonatal			•	•
989803142711	n/a	View >	Verification gas	1 gas cylinder	Mainstream CO ₂ sensor: M2501A, MP2: 865040; X2: 865039; X3: 867030, Xper Flex: 860335, 860338; Capnography extension 867040;	•	•



989803142651
M2501A



989803142661
M2533A



989803142671
M2536A



989803206601



989803206591



989803142711



Cardiac output

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803104611	M1642A	Cardiac output cable for use with transpulmonary thermodilution Length: 2,9m (9.5 ft)	1 cable	CM12: 863303, 863366; CM120: 863302, 863321; CM150: 863304, 863322 Hemo extension: 867039, M3012A; Capno extension: 867040; Microstream™ extension: 867041; PM6000: 867311, 867313, 867315;	•	•
989803104621	M1643A	Cardiac output cable for use with pulse contour analysis and transpulmonary thermodilution Length: 4.8 m (15.7 ft)	1 cable		•	•
989803219991		FloTrac Pressure Cable	1 cable	FloTrac module 867409	•	



989803104611
M1642A



989803104621
M1643A



989803219991



Electrodes

Single-patient use ECG monitoring electrodes

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA	
Wet Gel								
989803101671	40493D	View >	Wet gel ECG electrode, adult, foam Shelf life: 6 months, diameter: 55 mm (2.13 in)	5 electrodes/pouch 15 pouches/box 4 boxes/case (300 electrodes/case)	Adult ECG monitoring	•	•	
989803101681	40493E	View >	Wet gel ECG electrode, adult, foam Shelf life: 6 months, diameter: 54 mm (2.13 in)	30 electrodes/pouch 10 pouches/box (300 electrodes/box)		•	•	
989803192551		View >	Wet gel ECG electrode, adult, clear tape	30 electrodes/pouch 10 pouches/box (300 electrodes/box)		•	•	
Solid Gel								
989803129091	M4612A	View >	Solid gel ECG electrode, adult, foam, teardrop shape Shelf life: 6 months, diameter: 50 mm (1.97 in)	5 electrodes/pouch 60 pouches/box (300 electrodes/box)		•	•	
989803129101	M4613A	View >	Solid gel ECG electrode, adult, foam, teardrop shape Shelf life: 6 months, diameter: 50 mm (1.97 in)	30 electrodes/pouch 10 pouches/box (300 electrodes/box)		•	•	
Radiolucent								
989803148821		View >	Solid gel ECG electrode, small adult, radiolucent foam Shelf life: 6 months, size: 35 x 40 mm (1.4 x 1.6 in)	30 electrodes/pouch 20 pouches/box (600 electrodes/box)		•	•	
989803105971	M2202A	View >	Radiolucent electrode, adult Shelf life: 6 months, diameter: 50 mm (1.97 in)	30 electrodes/pouch 10 pouches/box (300 electrodes/box)		•	•	



989803101671
40493D



989803101681
40493E




989803192551




989803129091
M4612A



989803129101
M4613A



989803148821




989803105971
M2202A




Electrodes

Adult single-patient use pre-attached lead wire electrodes- AAMI color coding

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
Metallic						
989803156201		Pre-attached metallic 3-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in)	3 electrodes/set 1 set/pouch 100 pouches/box (300 electrodes/box)	Adult ECG monitoring	•	•
989803156241		Pre-attached metallic 5-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in)	5 electrodes/set 1 set/pouch 60 pouches/box (300 electrodes/box)		•	•
Radiolucent						
989803156221		Pre-attached radiolucent 3-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in)	3 electrodes/set 1 set/pouch 100 pouches/box (300 electrodes/box)		•	•
989803156261		Pre-attached radiolucent 5-lead wire (0.91 m / 36 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in)	5 electrodes/set 1 set/pouch 60 pouches/box (300 electrodes/box)		•	•



989803156201



989803156241

Electrodes

Infant/Pediatric single-patient use pre-attached lead wire electrodes- AAMI color coding


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
Metallic						
989803100501	13952A	Pre-attached metallic 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in); socket size 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)	Infant/pediatric ECG monitoring		•
989803100541	13953A	Pre-attached metallic 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)			•
Radiolucent						
989803100511	13952B	Pre-attached radiolucent 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in); socket size: 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)			•
989803100551	13953B	Pre-attached radiolucent 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size: 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)			•




989803100501
13952A



989803100541
13953A



989803100511
13952B




989803100551
13953B

Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.


Electrodes

NeoLead Infant/Pediatric single-patient use pre-attached lead wire electrodes- AAMI color coding (High humidity environment, 5-7days)


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA	
Metallic							
989803183101		NeoLead, metallic pre-attached 61 cm (24 in) 3-lead wire pediatric/neonatal ECG electrodes, hydrogel/hydrocolloid adhesive, AAMI, standard DIN connector Diameter: 25 mm (1 in)	3 sets/pouch 50 pouches/box (150 sets/box)	Infant/pediatric ECG monitoring		•	
989803183081		Mini Neolead, metallic pre-attached 61 cm (24 in) 3-lead wire pediatric/neonatal ECG electrodes, hydrogel/hydrocolloid adhesive, AAMI, standard DIN connector Diameter: 19 mm (0.75 in)	3 sets/pouch 50 pouches/box (150 sets/box)			•	
Radiolucent							
989803183121		NeoLead, radiolucent pre-attached 61 cm (24 in) 3-lead wire pediatric/neonatal ECG electrodes, AAMI, standard DIN connector Diameter: 25 mm (1 in)	3 sets/pouch 50 pouches/box (150 sets/box)			•	
989803183141		Micro NeoLead, radiolucent pre-attached 61 cm (24 in) 3-lead wire pediatric/neonatal ECG electrodes, AAMI, standard DIN connector Dimensions: 25 x 9.5 mm (1 x 0.37 in)	3 sets/pouch 50 pouches/box (150 sets/box)			•	




989803183101



989803183081



989803183121



989803183141

Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.



Electrodes

Adult single-patient use pre-attached lead wire electrodes- IEC color coding

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
Radiolucent						
989803156231		Pre-attached radiolucent 3-lead wire (1.0 m / 39.4 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in)	3 electrodes/set 1 set/pouch 100 pouches/box (300 electrodes/box)	Adult ECG monitoring: IntelliVue	•	
989803156271		Pre-attached radiolucent 5-lead wire (1.0 m / 39.4 in) electrode set, non-shielded, solid gel Shelf life: 24 months unopened, dimensions: 30 x 45 mm (1.18 x 1.77 in)	5 electrodes/set 1 set/pouch 60 pouches/box (300 electrodes/box)		•	



989803156231




989803156271

Electrodes

Infant/Pediatric single-patient use pre-attached lead wire electrodes- IEC color coding


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
Metallic						
989803100521	13952D	Pre-attached metallic 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in); socket size: 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)	infant/pediatric ECG monitoring	•	•
989803100561	13953D	Pre-attached metallic 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size: 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)		•	•
Radiolucent						
989803100571	13953E	Pre-attached radiolucent 46 cm (18 in) 3-lead wire electrode, non-shielded, solid gel Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in); socket size: 1.5 mm (0.06 in)	3 sets/pouch 100 pouches/box (300 sets/box)		•	•



989803100521
13952D



989803100561
13953D



989803100571
13953E

Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.



Cables and leads

3-lead monitoring - AAMI color coding

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803145091	M1671A	View >	3-lead Grabber, limb leads, reusable, length: 1 m (3.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803145111	M1673A	View >	3-lead Snap, limb leads, reusable, length 1 m (3.3 ft)	1 each		•	•
989803144941	M1624A	View >	3L Mini-clip, Unshielded, Grabber, reusable, length 0.7 m (2.3 ft)	1 each		•	•
989803145131	M1675A	View >	3-lead OR Grabber, limb leads, reusable, for use in operating room, length: 1 m (3.3 ft)	1 each		•	•
989803173121		View >	3-lead Grabber, limb leads, single patient use, length: 1 m (3.3 ft)	20 each		•	•
989803145071	M1669A	View >	3-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included 3 leads, reusable, length: 2.7 m (9.9 ft)	1 cable		•	•
989803170171		View >	3-lead OR trunk cable, orange, 12-pin connector, AAMI and IEC labels included For use in operating room 3 leads, reusable, length: 2.7 m (9.9 ft)	1 cable		•	•



989803145091
M1671A



989803145111
M1673A




989803144941
M1624A



989803145131
M1675A



989803173121



989803145071
M1669A




989803170171




Cables and leads

5-lead monitoring - AAMI color coding


Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803125841	M1968A	View >	5-lead Grabber, limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803144991	M1644A	View >	5-lead Snap, limb leads, reusable, Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each		•	•
989803125861	M1973A	View >	5-lead OR, grabber, limb leads, reusable, for use in operating room, Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each		•	•
989803173131		View >	5-lead Grabber, limb leads, single patient use, Length: 1 m (3.3 ft)	20 each		•	•
989803173151		View >	5-lead grabber w/ shower shield, limb leads, single-patient use, Length: 85 cm (34 in.)	20 each	TRx4851A: 862231; TRx4841A: 862439	•	•
989803145061	M1668A	View >	5-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included, 5 leads, reusable, Length: 2.7 m (9.9 ft)	1 cable	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803170181		View >	5-lead OR ECG trunk cable, orange, 12-pin connector, AAMI and IEC labels included. For use in operating room, 5 leads, reusable, Length: 2.7 m (9.9 ft)	1 cable		•	•



989803125841
M1968A



989803144991
M1644A



989803125861
M1973A



989803173131



989803173151



989803145061
M1668A




989803170181



Cables and leads

6-lead monitoring - AAMI color coding

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803145161	M1680A	6-lead Grabber, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803197511		6-lead Grabber, single-patient use Length: 1 m (3.3 ft)	20 each		•	•
989803145051	M1667A	6-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included, reusable Length: 2.7 m (9.9 ft)	1 cable		•	•



989803145161
M1680A



989803197511




989803145051
M1667A




Cables and leads

10-lead monitoring - AAMI color coding. (requires 5-lead limb + 5-lead chest)


Ordering code	Alternative ordering code		Description	Packaging	Use with	CE	FDA
limb leads							
989803125841	M1968A	View >	5-lead Grabber, limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803144991	M1644A	View >	5-lead Snap, limb leads, reusable, Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each		•	•
989803125861	M1973A	View >	5-lead OR, grabber, limb leads, reusable, for use in operating room, Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each		•	•
989803173131		View >	5-lead Grabber, limb leads, single patient use, Length: 1 m (3.3 ft)	20 each		•	•
989803173151		View >	5-lead grabber w/ shower shield, limb leads, single-patient use, Length: 85 cm (34 in.)	20 each	TRx4851A: 862231; TRx4841A: 862439	•	•
chest leads							
989803125881	M1976A	View >	5-lead Grabber, chest leads, reusable, Length: 1 m (3.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 5+5 trunk cable M1949A, M1663A	•	•
989803144911	M1602A	View >	5-lead Snap, chest leads, reusable, Length: 1 m (3.3 ft)	1 each		•	•
989803125901	M1979A	View >	5-lead OR Grabber, chest leads, reusable, for use in operating room, Length: 1 m (3.3 ft)	1 each		•	•
989803192141		View >	5-lead grabber, chest leads, single-patient use, Length: 85 cm (34 in.)	20 each		•	•
cables							
989803125831	M1949A	View >	10-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included. 10 (5-lead limb + 5-lead chest) leads, reusable, Length: 2.7 m (9.9 ft)	1 cable		•	•
989803144791	M1663A	View >	10-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included. 10 (5-lead limb + 5-lead chest) leads, reusable, Length: 2 m (6.5 ft)	1 cable		•	•



989803125841
M1968A



989803144991
M1644A




989803125861
M1973A



989803173131



989803173151



989803125881
M1976A



989803144911
M1602A



989803125901
M1979A



989803192141



989803125831
M1949A



989803144791
M1663A



Cables and leads

10-lead monitoring - AAMI color coding (requires 6-lead limb + 4-lead chest).

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
limb leads						
989803145161	M1680A	6-lead Grabber, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803197511		6-lead Grabber, single-patient use Length: 1 m (3.3 ft)	20 each		•	•
chest leads						
989803144841	M1532A	4-lead Grabber, chest leads, reusable, length: 1 m (3.3ft) Must be ordered in pairs with 989803145161 (M1680A)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 6+4 trunk cable M1665A, PM6000: 867311, 867313, 867315	•	•
cable						
989803145041	M1665A	10-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included. 10 (6-lead limb + 4-lead chest) leads, reusable. Length: 2.7 m (9.9 ft)	1 cable	MP2: 865040; MP5: 865024; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033, PM6000: 867311, 867313, 867315	•	•



989803145161
M1680A



989803197511



989803144841
M1532A




989803145041
M1665A



Cables and leads

3-lead monitoring - IEC color coding


Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803145101	M1672A	View >	3-lead Grabber, limb leads, reusable Length: 1 m (3.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803145121	M1674A	View >	3-lead Snap, limb leads, reusable Length: 1 m (3.3 ft)	1 each		•	•
989803144951	M1626A	View >	3-lead Miniclip, Unshielded, Grabber, reusable, length 0.7 m (2.3 ft)			•	•
989803145141	M1678A	View >	3-lead OR, Grabber, limb leads, reusable, for use in operating room Length: 1 m (3.3 ft)	1 each		•	•
989803174201		View >	3-lead Grabber, limb leads, single patient use Length: 1 m (3.3 ft)	20 each		•	•
989803145071	M1669A	View >	ECG 3-lead trunk cable, 12-pin connector, AAMI and IEC labels included 3 leads, reusable, length: 2.7 m (9.9 ft)	1 cable	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; 3-ch ECG module: 862275; 3-ch ECG/respiration module: 862276; VM3: 863077; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; DFM100: 866199; HeartStart XL: 861266; HeartStart XL+: 861290; HeartStart MRx: 861288, 861289	•	•
989803170171		View >	3-lead OR ECG trunk cable, orange, 12-pin connector, AAMI and IEC labels included For use in operating room 3 leads, reusable, length: 2.7 m (9.9 ft)	1 cable	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; VM6: 863064, 863065; VM8: 863066, 863068, 863317; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100: 866199	•	•



989803145101
M1672A



989803145121
M1674A



989803144951
M1626A



989803145141
M1678A



989803174201



989803145071
M1669A




989803170171




Cables and leads

5-lead monitoring - IEC color coding


Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803125851	M1971A	View >	5-lead grabber, limb leads, reusable Length: limb 1.6 m (5.3 ft), chest 1 m (3.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803145001	M1645A	View >	5-lead snap, limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each		•	•
989803125871	M1974A	View >	5-lead OR grabber, limb leads, reusable, for use in operating room Length: limb 1.6 m (5.3 ft), chest 1 m (3.3 ft)	1 each		•	•
989803174211		View >	5-lead grabber, limb leads, single-patient use Length: 1 m (3.3 ft)	20 each		•	•
989803145061	M1668A	View >	5-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included 5 leads, reusable Length: 2.7 m (9.9 ft)	1 cable		•	•
989803170181		View >	5-lead OR ECG trunk cable, orange, 12-pin connector, AAMI and IEC labels included For use in operating room 5 leads, reusable Length: 2.7 m (9.9 ft)	1 cable		•	•



989803125851
M1971A



989803145001
M1645A



989803125871
M1974A



989803174211



989803145061
M1668A




989803170181



Cables and leads

6-lead monitoring - IEC color coding

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803145171	M1681A	6-lead grabber, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803197521		6-lead grabber, single patient use Length: 1 m (3.3 ft)	20 each		•	•
989803145051	M1667A	6-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included 6 leads, reusable Length: 2.7 m (9.9 ft)	1 cable	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MX100: 867033; HeartStart XL+: 861290, PM6000: 867311, 867313, 867315	•	•



989803145171
M1681A



989803197521




989803145051
M1667A




Cables and leads

10-lead monitoring - IEC color coding. (requires 5-lead limb+ 5-lead chest)


Ordering code	Alternative ordering code		Description	Packaging	Use with	CE	FDA
limb leads							
989803125851	M1971A	View >	5-lead grabber, limb leads, reusable Length: limb 1.6 m (5.3 ft), chest 1 m (3.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803145001	M1645A	View >	5-lead snap, limb leads, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each		•	•
989803125871	M1974A	View >	5-lead OR grabber, limb leads, reusable, for use in operating room Length: limb 1.6 m (5.3 ft), chest 1 m (3.3 ft)	1 each		•	•
		View >	5-lead grabber, limb leads, single-patient use Length: 1 m (3.3 ft)	20 each		•	•
chest leads							
989803125891	M1978A	View >	5-lead grabber, chest leads, reusable Length: 1 m (3.3 ft)	1 each		•	•
		View >	5-lead grabber, chest leads, single patient use Length: 1 m (3.3 ft)	20 each		•	•
Cables							
989803125831	M1949A	View >	10-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included. 10 (5-lead limb + 5-lead chest) leads, reusable. Length: 2.7 m (9.9 ft)	1 cable		•	•
989803144791	M1663A	View >	10-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included. 10 (5-lead limb + 5-lead chest) leads, reusable. Length: 2 m (6.5 ft)	1 cable		•	•



989803125851
M1971A




989803145001
M1645A




989803125871
M1974A



989803174211



989803125891
M1978A



989803192151



989803125831
M1949A




989803144791
M1663A



Cables and leads

10-lead monitoring - IEC color coding (requires 6-lead limb + 4-lead chest)

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
limb leads						
989803145171	M1681A	6-lead grabber, reusable Length: limb 1 m (3.3 ft), chest 1.6 m (5.3 ft)	1 each	MP2: 865040; MP5: 865024; MMS: 862392, 862442; X2: 865039; X3: 867030; MMX: 867033; PM6000: 867311, 867313, 867315; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; DFM100:866199;	•	•
989803197521		6-lead grabber, single patient use Length: 1 m (3.3 ft)	20 each		•	•
chest leads						
989803144851	M1533A	4-lead Grabber, chest leads, reusable Length: 1 m (3.3 ft). Must be ordered in pairs with 989803145171 (M1681A)	1 each		•	•
cable						
989803145041	M1665A	10-lead ECG trunk cable, 12-pin connector, AAMI and IEC labels included. 10 (6-lead limb + 4-lead chest) leads, reusable. Length: 2.7 m (9.9 ft)	1 cable	MP2: 865040; MP5: 865024; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033, PM6000: 867311, 867313, 867315	•	•



989803145171
M1681A



989803197521



989803144851
M1533A



989803145041
M1665A



Cables and leads

Efficia

3-lead sets

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803160651		Grabber, limb leads, AAMI, reusable	1 each	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322, PM6000: 867311, 867313, 867315	•	
989803160661		Grabber, limb leads, IEC, reusable	1 each		•	

5-lead sets

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803160691		Grabber, limb leads, AAMI, reusable	1 each	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322, PM6000: 867311, 867313, 867315	•	
989803160701		Grabber, limb leads, IEC; reusable	1 each		•	

3/5-lead sets

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803160641		3/5 leads ECG trunk cable, AAMI/IEC, reusable	1 each	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322, PM6000: 867311, 867313, 867315	•	




989803160651



989803160661



989803160691



989803160701



989803160641



Cables and Leads

Telemetry Lead Sets and Cable – AAMI color coding

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803171811		View >	3-lead ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in)	1 each	MX40: 865350, 865351, 865352	•	•
989803171821		View >	5-lead ECG, snaps, reusable, length: 85 cm (34 in)	1 each		•	•
989803171831		View >	5-lead ECG, grabbers, reusable, length: 85 cm (34 in)	1 each		•	•
989803171841		View >	5-lead ECG + SpO ₂ , snaps, reusable, length: 85 cm (34 in)	1 each		•	•
989803171851		View >	5-lead ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in)	1 each		•	•
989803171871		View >	6-lead ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in)	1 each		•	•
989803172031		View >	5-lead ECG, grabbers, single-patient use, length 85 cm (34 in)	20 each		•	•
989803172051		View >	5-lead ECG + SpO ₂ , grabbers, single-patient use, length: 85 cm (34 in) FAST Only	20 each		•	•
989803202411		View >	5-lead ECG, grabber AAMI, Masimo SpO ₂ , telemetry Length: 85 cm (34 in)	1 each		•	•
989803172221		View >	MX40 Telemetry ECG trunk cable, AAMI/IEC Length: 2.7 m (9.9 ft)			•	•



989803171811



989803171821



989803171831



989803171841



989803171851



989803171871



989803172031



989803172051



989803202411



989803172221




Cables and Leads

Telemetry Lead Sets and Cable – IEC color coding


Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803171911			3-lead ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in)	1 each	MX40: 865350, 865351, 865352, M6000: 867311, 867313, 867315	•	•
989803171931			5-lead ECG, grabbers, reusable, length: 85 cm (34 in)	1 each		•	•
989803171951			5-lead ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in)	1 each		•	•
989803171971			6-lead ECG + SpO ₂ , grabbers, reusable, length: 85 cm (34 in)	1 each		•	•
989803172131			ECG, 5 leads, grabbers, single-patient use, length: 85 cm (34 in)	20 each		•	•
989803172151		View >	5-lead ECG + SpO ₂ , grabbers, single-patient use, length: 85 cm (34 in)	20 each		•	•
989803202421			5-lead ECG, grabber, IEC, Masimo SpO ₂ , telemetry length: 85 cm (34 in)	1 each		•	•
989803172221		View >	MX40 Telemetry ECG trunk cable, AAMI/IEC, Length: 2.7 m (9.9 ft)	1 each	MX40: 865350, 865351, 865352	•	•




989803171911



989803171931



989803171951



989803171971



989803172131



989803172151



989803172221



Telemetry Adapters

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
Telemetry FAST SpO ₂							
989803199081		View >	MX40 telemetry adapter, SpO ₂ only	1 each	MX40: 865350, 865351, 865352	•	•
989803199101		View >	MX40 telemetry adapter, ECG + SpO ₂ , 3 leads	1 each		•	•
989803199091		View >	MX40 telemetry adapter, ECG + SpO ₂ , 5 leads	1 each		•	•
989803199071		View >	MX40 telemetry adapter, ECG + SpO ₂ , 6 leads	1 each		•	•
Telemetry Masimo SpO ₂							
989803202391		View >	MX40 Adapter Masimo SpO ₂ -only	1 each	MX40: 865350, 865351, 865352	•	•
989803202381		View >	MX40 Adapter 5-lead + Masimo SpO ₂	1 each		•	•
989803204291		View >	MX40 Adapter 6-lead + Masimo SpO ₂	1 each		•	•



989803199081



989803199101



989803199091



989803199071



989803202391



989803202381



989803204291



Trunk Cable Adapters (Non-Philips monitors)

Space Lab Ultra ECG Adapters

Set	Philips Reorder Code	Image	Description	Use with	Reusable/disposable	Monitor connectivity	CE	FDA
3-lead	989803211121	View >	Spacelabs Ultra-Philips		Reusable	UltraView SL		•
5-lead	989803211131	View >	Spacelabs Ultra-Philips		Reusable	UltraView SL		•

Mindray Pano ECG Adapters

Set	Philips Reorder Code	Image	Description	Use with	Reusable/disposable	Monitor connectivity	CE	FDA
3-lead	989803211141	View >	Mindray Pano-Philips		Reusable	Panorama		•
5-lead	989803211151	View >	Mindray Bene/Pass-Philips		Reusable	Beneview/Passport		•

Nihon Kohden ECG Adapters

Set	Philips Reorder Code	Image	Description	Use with	Reusable/disposable	Monitor connectivity	CE	FDA
3-lead	989803211161	View >	Nihon Kohden-Philips		Reusable	Lifescopes/BSM Series		•
6-lead	989803211171	View >	Nihon Kohden-Philips		Reusable	Lifescopes/BSM Series		•

GE-Philips ECG Adapters

Set	Philips Reorder Code	Image	Description	Use with	Reusable/disposable	Monitor connectivity	CE	FDA
3-lead	989803211081	View >	GE-Philips		Reusable	CareScape B Series		•
5-lead	989803211091	View >	GE-Philips		Reusable	CareScape B Series / MAC 3500, 2000, 800, 600 ECG carts		•
5-lead	989803211111	View >	GE Apex Pro FH-Philips		Reusable	Apex Pro FH Telemetry		•
12-lead	989803211101	View >	GE-Philips		Reusable	CareScape B Series		•



989803211121



989803211131



989803211141



989803211151



989803211161



989803211171



989803211081



989803211091



989803211111



989803211101



Electrodes

Single-patient use pediatric EEG monitoring electrodes

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803100491	13951C	Disposable, repositionable solid gel electrode, snap 3 electrodes mounted on plastic card Shelf life: 6 months, dimensions: 22 x 33 mm (0.87 x 1.30 in)	3 electrodes/pouch 100 pouches/box (300 electrodes/box)	EEG module: 866078	•	•
989803100581	13955C	Disposable, repositionable solid gel electrode, snap 3 electrodes mounted on plastic card Shelf life: 6 months, dimensions: 22 x 22 mm (0.87 x 0.87 in)	3 electrodes/pouch 100 pouches/box (300 electrodes/box)		•	•




989803100491
13951C

989803100581
13955C

Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.

EEG supplies

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803105621	M1931A		5-lead EEG adult cup electrode set, reusable, length: 80 cm (31.5 in)	5 leads	EEG module: 866078	•	•
989803105631	M1932A		5-lead EEG pediatric cup electrode set, reusable, length 80 cm (31.5 in)	5 leads		•	•
989803105641	M1934A		5-lead EEG miniclip lead set, reusable, length: 80 cm (31.5 in)	5 leads		•	•
989803106091	M2268A		5-lead trunk cable, reusable, length: 2.7 m (8.8 ft)	1 cable		•	•
989803180521			9-lead EEG infant and pediatric cup electrode set, reusable, length: 80 cm (31.5 in)	9 leads		•	•
989803180531		View >	9-lead EEG miniclip lead set, reusable, length: 80 cm (31.5 in)	9 leads		•	•
989803180541		View >	4-channel trunk cable, reusable, length: 2.7 m (8.8 ft)	1 cable		•	•
989803180561		View >	9-lead EEG adult cup electrode set, reusable, length: 80 cm (31.5 in)	9 leads		•	•



989803105621
M1931A



989803105631
M1932A




989803105641
M1934A



989803106091
M2268A



989803180521



989803180531



989803180541



989803180561



IntelliVue NMT

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803174581	n/a	Patient cable with 12-pin instrument connector, reusable, length: 3.0 m (9.8 ft)	1 each	NMT module: 865383	•	•
989803199211	n/a	IntelliVue NMT hand adapter	10 pieces/box		•	•



989803174581



989803199211



Reusable cuffs

Easy care multi-patient cuffs - single hose

Philips reorder code	Size	Image	Description	Color	CE	FDA
989803147841 M4552B	Infant		Bladder length: 15 cm (5.9 in), bladder width: 5.5 cm (2.2 in), limb circumference: 10-15 cm (3.9-5.9 in) Packaging 1 cuff	Orange	•	•
989803147851 M4553B	Pediatric		Bladder length: 21.5 cm (8.5 in), bladder width: 8 cm (3.1 in), limb circumference: 14-21.5 cm (5.5-8.5 in) Packaging 1 cuff	Green	•	•
989803147861 M4554B	Small adult		Bladder length: 24 cm (9.4 in), bladder width: 10.6 cm (4.2 in), limb circumference: 20.5-28.5 cm (8-11.2 in) Packaging 1 cuff	Light blue	•	•
989803147871 M4555B	Adult		Bladder length: 30 cm (11.8 in), bladder width: 13.5 cm (5.3 in), limb circumference: 27.5-36 cm (10.8-14.1 in) Packaging 1 cuff	Navy blue	•	•
989803147881 M4556B	Adult extra-long*		Bladder length: 30 cm (11.8 in), bladder width: 13.5 cm (5.3 in), limb circumference: 27.5-36 cm (10.8-14.1 in) Packaging 1 cuff	Navy blue	•	•
989803147891 M4557B	Large adult		Bladder length: 38 cm (15 in), bladder width: 17 cm (6.7 in), limb circumference: 35-45 cm (13.7-17.7 in) Packaging 1 cuff	Burgundy	•	•
989803147901 M4558B	Large adult extra-long**		Bladder length: 38 cm (15 in), bladder width: 17 cm (6.7 in), limb circumference: 35-45 cm (13.7-17.7 in) Packagin 1 cuff	Burgundy	•	•
989803147911 M4559B	Thigh		Bladder length: 48 cm (18.9 in), bladder width: 21 cm (8.3 in), limb circumference: 44-56 cm (17.3-22 in) Packaging 1 cuff	Grey	•	•



Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.



Reusable cuffs

Comfort care multi-patient cuffs - single hose

Philips reorder code	Size	Image	Description	Color	CE	FDA
989803104141 M1571A	Infant		Bladder length: 15 cm (5.9 in), bladder width: 5.5 cm (2.2 in), limb circumference: 10-15 cm (3.9-5.9 in) Packaging 1 cuff	Brown	•	•
989803104151 M1572A	Pediatric		Bladder length: 21.5 cm (8.5 in), bladder width: 8 cm (3.1 in), limb circumference: 14-21.5 cm (5.5-8.5 in) Packaging 1 cuff	Brown	•	•
989803104161 M1573A	Small adult		Bladder length: 24 cm (9.4 in), bladder width: 10.5 cm (4.1 in), limb circumference: 20.5-28 cm (8-11 in) Packaging 1 cuff	Brown	•	•
989803104171 M1574A	Adult		Bladder length: 28 cm (11 in), bladder width: 13 cm (5.1 in), limb circumference: 27-35 cm (10.6-13.7 in) Packaging 1 cuff	Brown	•	•
989803191311 M1574XL	Adult extra-long*		Bladder length: 28 cm (11 in), bladder width: 13 cm (5.1 in), limb circumference: 27-35 cm (10.6-13.7 in) Packaging 1 cuff	Brown	•	•
989803104181 M1575A	Large adult		Bladder length: 35 cm (13.8 in), bladder width: 16 cm (6.3 in), limb circumference: 34-43 cm (13.3-16.9 in) Packaging 1 cuff	Brown	•	•
989803104191 M1576A	Thigh		Bladder length: 43 cm (16.9 in), bladder width: 20 cm (7.9 in), limb circumference: 42-54 cm (16.5-21.2 in) Packaging 1 cuff	Brown	•	•

Accuracy of extra-long cuffs is dependent on proper application according to cuff labeling.
*26.0 inches (66.0 cm) total cuff length.
**32.9 inches (83.5 cm) total cuff length"



Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.



Disposable cuffs

Gentle care single patient disposable NiBP cuffs

Philips reorder code	Size	Image	Description	Color	CE	FDA
989803148001 M4572B	Infant		Bladder length: 15 cm (5.9 in), bladder width: 5.5 cm (2.2 in), limb circumference: 10-15 cm (3.9-5.9 in) Packaging 10 cuffs	Orange	•	•
989803148011 M4573B	Pediatric		Bladder length: 21.5 cm (8.5 in), bladder width: 8 cm (3.1 in), limb circumference: 14-21.5 cm (5.5-8.5 in) Packaging 10 cuffs	Green	•	•
989803148021 M4574B	Small adult		Bladder length: 24 cm (9.4 in), bladder width: 10.6 cm (4.2 in), limb circumference: 20.5-28.5 cm (8-11.2 in) Packaging 10 cuffs	Light blue	•	•
989803148031 M4575B	Adult		Bladder length: 30 cm (11.8 in), bladder width: 13.5 cm (5.3 in), limb circumference: 27.5-36 cm (10.8-14.1 in) Packaging 10 cuffs	Navy blue	•	•
989803148041 M4576B	Adult extra-long*		Bladder length: 30 cm (11.8 in), bladder width: 13.5 cm (5.3 in), limb circumference: 27.5-36 cm (10.8-14.1 in) Packaging 10 cuffs	White/navy with striping	•	•
989803148051 M4577B	Large adult		Bladder length: 38 cm (15 in), bladder width: 17 cm (6.7 in), limb circumference: 35-45 cm (13.7-17.7 in) Packaging 10 cuffs	Burgundy	•	•
989803148061 M4578B	Large adult extra-long**		Bladder length: 38 cm (15 in), bladder width: 17 cm (6.7 in), limb circumference: 35-45 cm (13.7-17.7 in) Packaging 10 cuffs	Burgundy with striping	•	•




Accuracy of extra-long cuffs is dependent on proper application according to cuff labeling.
*26.0 inches (66.0 cm) total cuff length.
**32.9 inches (83.5 cm) total cuff length”



Neonatal products

Soft neonatal/infant cuffs - single hose

Philips reorder code	Size	Image	Description	Noenatal/infant	Color	CE	FDA
989803167231 M1866S	#1		Limb circumference: 3.1-5.7 cm (1.2-2.2 in) Packaging 40 cuffs	Neonatal	Orange	•	•
989803167241 M1868S	#2		Limb circumference: 4.3-8.0 cm (1.7-3.1 in) Packaging 40 cuffs	Neonatal	Light blue	•	•
989803167251 M1870S	#3		Limb circumference: 5.8-10.9 cm (2.2-4.3 in) Packaging 40 cuffs	Neonatal	Green	•	•
989803167261 M1872S	#4		Limb circumference: 7.1-13.1 cm (2.8-5.5 in) Packaging 40 cuffs	Neonatal	Navy blue	•	•
989803167271 M1873S	#5		Limb circumference: 10-15 cm (3.9-5.9 in) Packaging 40 cuffs	Infant	Burgundy	•	•
989803167551		View >	Single-patient NiBP Kit Packaging 50 cuffs, consists of: 5x M1866S, 10x M1868S, 20x M1870S, 10x M1872S, 5x M1873S	Neonatal	Various	•	•



989803167231
M1866S



989803167241
M1868S



989803167251
M1870S



989803167261
M1872S



989803167271
M1873S



989803167551



Neonatal products

Vinyl neonatal/infant cuffs - single hose

Philips reorder code	Size	Description	Noenatal/infant	Color	CE	FDA
989803167181 M1866B	#1	Limb circumference: 3.1-5.7 cm (1.2-2.2 in) Packaging 40 cuffs	Neonatal	Orange	•	•
989803167191 M1868B	#2	Limb circumference: 4.3-8.0 cm (1.7-3.1 in) Packaging 40 cuffs	Neonatal	Light blue	•	•
989803167201 M1870B	#3	Limb circumference: 5.8-10.9 cm (2.2-4.3 in) Packaging 40 cuffs	Neonatal	Green	•	•
989803167211 M1872B	#4	Limb circumference: 7.1-13.1 cm (2.8-5.5 in) Packaging 40 cuffs	Neonatal	Navy blue	•	•
989803167221 M1873B	#5	Limb circumference: 10-15 cm (3.9-5.9 in) Packaging 40 cuffs	Infant	Burgundy	•	•
989803167541		Single-patient NiBP Kit Packaging 50 cuffs, consists of: 5x M1866S, 10x M1868S, 20x M1870S, 10x M1872S, 5x M1873S	Neonatal	Various	•	•





Accessories and hoses

New NiBP Hoses (use with 5mm NiBP monitor connection)

Philips reorder code	Description	CE	FDA
989803209761	Adult NiBP air hose Length: 1.5 m (4.92 ft) Packaging 1 air hose	•	•
989803209771	Adult NiBP air hose Length: 3.0 m (9.8 ft) Packaging 1 air hose	•	•
989803209741	Neonate NiBP air hose Length: 1.5 m (4.92 ft) Packaging 1 air hose	•	•
989803209751	Neonate NiBP air hose Length: 3.0 m (9.8 ft) Packaging 1 air hose	•	•

Ask your Philips representative about the need for Conversion Kits when changing from Traditional Hoses to New Hoses



989803209761



989803209771



989803209741



989803209751



Accessories and hoses

Traditional NiBP Hoses (use with 6mm NiBP monitor connection)

Philips reorder code	Description	CE	FDA
989803104331 M1598B	Adult NiBP air hose Length: 1.5 m (4.92 ft) Packaging 1 air hose	•	•
989803104341 M1599B	Adult NiBP air hose Length: 3.0 m (9.8 ft) Packaging 1 air hose	•	•
989803166851 M1596C	Neonate NiBP air hose Length: 1.5 m (4.92 ft) Packaging 1 air hose	•	•
989803166861 M1597C	Neonate NiBP air hose Length: 3.0 m (9.8 ft) Packaging 1 air hose	•	•

Ask your Philips representative about the need for Conversion Kits when changing from Traditional Hoses to New Hoses



989803104331
M1598B



989803104341
M1599B



989803166851
M1596C



989803166861
M1597C



Conversion kit

New ISO 80369-1 Compliant NBP Inlet Conversion Kit

Philips reorder code	Description	Packaging	Use with	CE	FDA
989803217001	NBP inlet 5 mm	5 each	Compatible with all monitors except X3, MX100, MMX	•	•
989803217011	NBP inlet 5 mm	50 each	Compatible with all monitors except X3, MX100, MMX	•	•
989803217021	NBP inlet covered 5 mm	5 each	Compatible with X3, MX100, MMX	•	•
989803217031	NBP inlet covered 5 mm	50 each	Compatible with X3, MX100, MMX	•	•



Reusable invasive blood pressure

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803102481	CPJ84022	View >	Single-use sterile domes	50 domes	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304; X3: 867030; MX100: 867033; MMX: 867036; Hemodynamic ext.: 862111, 867039; Capnography ext.: 862187, 867040; Microstream™ ext.: 867041 862393, 865377; X2: 865039; MP2: 865040; PM6000: 867313, 867315; IBP module 862277, MMS: 862442; MP5:865024;	•	•
989803102491	CPJ84046	View >	Transducer holder	4 holders		•	•
989803102511	CPJ840J6	View >	Reusable pressure transducer, used to measure arterial and venous blood pressure Includes one transducer holder Cable length: 3 m (9.8 ft)	1 transducer		•	•
989803102521	CPJ84447	View >	IV pole mount	1 mount		•	•
989803136281	CPJ84026	View >	Monitoring line set	10 sets/bag		•	•
989803199741		View >	Dual IBP adapter	1 adapter		•	•
989803220201		View >	TruWave Reusable Cable	1 cable		•	




989803102481



989803102491



989803102511



989803102521



989803136281



989803199741



989803220201



Recording paper

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803109871	M4816A	Thermal array recorder paper, 50 mm x 20 m (2 in x 66 ft)	10 rolls	Thermal Array Recorder Module M1116C, PM6000: 867311, 867313, 867315	•	•
989803136891	n/a	Thermal array recorder paper, 5.8 cm x 30.5 m (2.2 in x 100 ft)	5 rolls	CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304 SureSigns: 863277, 863278, 863279, 863283; VS30: 863359, 863380	•	•
989803190381	n/a	Efficia chemical thermal recorder paper, 50 mm (1.9 in), no lines	12 rolls	DFM 100: 866199	•	•



989803109871



989803136891



989803190381



Reusable sensors

Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803144371	M1191B	View >	Reusable adult SpO ₂ sensor, 8-pin connector Length: 2 m (6.5 ft) Patient size: > 50 kg (> 110 lb)	1 sensor	Pulse oximetry module: 862112; MMS: 862442, 866462; MP2: 865040; MP5: 865024, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 867146; MX100: 867033; MMX: 867036; VS2: 863278, 863279; VS4: 863283; VSi: 863276; 863380 CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; PM6000: 867311, 867313, 867315; FM30: 862199; FM40: 865132; FM50: 865071. Except product M1191T, M1196T: DFM100: 866199	•	•
989803144381	M1191BL	View >	Reusable adult SpO ₂ sensor, 8-pin connector Length: 3 m (9.8 ft) Patient size: > 50 kg (> 110 lb)	1 sensor		•	•
989803128591	M1191T	View >	Reusable adult SpO ₂ sensor, 9-pin (D-sub) connector Length: 0.45 m (1.5 ft) Patient size: > 50 kg (> 110 lb)	1 sensor		•	•
989803205871	M1192A	View >	Reusable pediatric/small adult finger SpO ₂ Sensor, 8-pin connector Length: 1.5 m (4.9 ft) Patient size: 15–50 kg (33–110 lb)	1 sensor		•	•
989803205881	M1193A	View >	Reusable neonatal SpO ₂ wrap sensor, 8-pin connector Length: 1.5 m (4.9 ft) Patient size: 1-4 kg	1 sensor		•	•
989803205891	M1194A	View >	Pediatric/adult ear clip SpO ₂ sensor, 8-pin connector Length: 1.5 m (4.9 ft) Patient size: > 40 kg (> 88 lb)	1 sensor		•	•
989803205901	M1196A	View >	Reusable adult SpO ₂ clip sensor, 8-pin connector Length: 3 m (9.8 ft) Patient size: > 40 kg (> 88 lb)	1 sensor		•	•
989803205911	M1196T	View >	Reusable adult SpO ₂ clip sensor, 9-pin (D-sub) connector Length: 0.90 m (2.9 ft) Patient size: >40 kg (> 88 lb)	1 sensor		•	•



989803144371
M1191B



989803144381
M1191BL



989803128591
M1191T



989803205871
M1192A




989803205881
M1193A



989803205891
M1194A



989803205901
M1196A



989803205911
M1196T



Single-patient sensors


Ordering code	Alternative ordering code	Image	Description	Packaging	Use with	CE	FDA
989803205381		View >	Nasal Alar SpO ₂ sensor, compatible with FAST technology	20 sensors	Pulse oximetry module: 862112; MMS: 862442; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 865352; MX100: 867033; MMX: 867036; MP2: 865040 FM30: 862199; FM40: 865132; FM50: 865071; PM6000: 867311, 867313, 867315	•	•
989803205391		View >	Nasal Alar SpO ₂ sensor, compatible with FAST and Nellcor™ technologies	24 sensors	Pulse oximetry module: 862112; MMS: 862442; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 865352; MX100: 867033; MMX: 867036; MP2: 865040 Nellcor N-600x; Nellcor N-595; DataScope Passport 2, FM30: 862199; FM40: 865132; FM50: 865071	•	•
989803205831	M1131A	View >	Disposable adult/pediatric SpO ₂ sensor, 9-pin (D-sub) connector Length: 45 cm (17.7 in) Patient size: > 20 kg (> 44 lb)	20 sensors	Pulse oximetry module: 862112; MMS: 862442, MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX40: 865350, 865351, 867146; MX100: 867033; MMX: 867036; VS2: 863278, 863279; VSi: 863276, 863277; VS30: 863380, CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322, FM30: 862199; FM40: 865132; FM50: 865071; XPER: 860338 Intrepid: 867172	•	•
989803205841	M1132A	View >	Infant disposable SpO ₂ sensor, 9-pin (D-sub) connector Length: 0.9 m (2.9 ft) Patient sizes: Infant (finger/toe): 3-10kg (6.6-22 lb)	20 sensors	Pulse oximetry module: 862112; MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; VM1: 863264, 863265, 863266; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS4: 863283; VSi: 863276; 863277; VS30: 863359, 863380; Intrepid: 867172; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322, M1132A, Except product M1134A: 860338; Except M1133A/M1134A: 860335; FM30: 862199; FM40: 865132; FM50: 865071; MX40: 865350, 865351, 867146.	•	•
989803205851	M1133A	View >	Neonatal/infant/adult disposable SpO ₂ sensor, 9-pin (D-sub) connector, Length: 0.90 m (2.9 ft) Patient sizes: Adult (finger): > 40 kg (> 88 lb), Neonate (foot/hand): < 3 kg (< 6.6 lb) Infant (big toe/thumb): 10–20 kg (22–44 lb)	20 sensors		•	•
989803205861	M1134A	View >	Adhesive-free neonatal/infant/adult disposable SpO ₂ sensor, 9-pin (D-sub) connector, Length: 0.90 m (2.9 ft) Patient sizes: Adult (finger): > 40 kg (> 88 lb), Neonate (foot/hand): < 3 kg (< 6.6 lb), Infant (big toe/thumb): 10–20 kg (22–44 lb)	20 sensors		•	•



989803205381
989803205391



989803205831
M1131A



989803205841
M1132A



989803205861
989803205851
M1133A, M1134A




Adapter cables

Adapter cables for Philips sensors


Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803105681	M1941A	SpO ₂ Extension Cable (8-pin), adds additional length between 8-pin sensors and 8-pin sockets Length: 2 m (6.6 ft)	1 cable	Pulse oximetry module: 862112; MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; VM4: 863063, 863085; VM6: 863064, 863065, 863086; VM8: 863066, 863068, 863087, 863088, 863317; VS2: 863278, 863279; VS4: 863283; VSi: 863276; VS30: 863359; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; HeartStart XL+: 861290 FM30: 862199; FM40: 865132; FM50: 865071 Except product M1941A, M1943A/AL: 863277, 865350, 865351, 867146, 867172, M2703A, M2704A, M2705A Except product M1941A only: 863310, 863311, 863276, 863359 Except product M1943A only: 860335, 860338 Except product M1943A/AL only: 863380Except product M1943A, M1943AL, M1143AL: DFM100: 866199 Except product M1943AL: HeartStart MRx: 861288, 861289, 861464, 861482, 861491, 860396, 860397, 860398, 861465	•	•
989803105691	M1943A	SpO ₂ Adapter Cable (8-pin), adapts 9-pin sensors to 8-pin sockets Length: 1.1 m (3.6 ft)	1 cable		•	•
989803128651	M1943AL	SpO ₂ Adapter Cable (8-pin), adapts 9-pin sensors to 8-pin sockets Length: 3 m (9.8 ft)	1 cable		•	•



989803105681



989803105691



989803128651



Adapter cables

Medtronic Nellcor Oximax™

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803136591	M1943NL	SpO ₂ Adapter Cable for use with Nellcor OxiMax™ SpO ₂ sensors, 8-pin to 9-pin D-sub adapter cable, 3 m (9.8 ft)	1 cable	Nellcor Oximax™ SpO ₂ sensors on Philips devices.	•	•



989803136591
M1943NL



Adapter cables

Masimo® rainbow® SET™ SpO₂ adapter cables

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803182411	n/a	rainbow® SET™ cable for rainbow, angled Length: 1.2 m (4 ft)	1 cable	Masimo® rainbow® Module: 867191; MMS: 866462; X2: 865039; X3: 86703; MMX: 867036; MP2: 865040; MP5: 865024, 865322	•	•
989803186471	n/a	rainbow® SET™ cable for rainbow, straight Length: 3.7 m (12 ft)	1 cable	Masimo® rainbow® Module: 867191; MMS: 866462; X2: 865039; X3: 86703; MMX: 867036; MP2: 865040; MP5 865024, 865322; VS4: 863283	•	•
989803186501	n/a	rainbow® SET™ cable for LNCS® Length: 3 m (10 ft)	1 cable		•	•
989803186531	n/a	rainbow® SET™ cable for rainbow, angled Length: 1.2 m (4 ft)	1 cable	MMS: 866462; MP2: 865040; MP5: 865024, 865322; X2: 865039; X3: 86703; MMX: 867036; Masimo® rainbow® Module: 867191; VS4: 863283; VS30: 863359, 863380; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150; 863304, 863322	•	•



989803182411



989803186471



989803186501



989803186531



Adapter cables

Masimo® rainbow® SET™ LNCS® SpO₂ adapter cables

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803186501	n/a	rainbow® SET™ cable for LNCS® Length: 3 m (10 ft)	1 cable	Masimo® rainbow® Module: 867191; MMS: 866462; X2: 865039; MMX: 867036; MP2: 865040; MP5: 865024, 865322; VS4: 863283	•	•
989803186531	n/a	rainbow® SET™ cable for LNCS®, straight Length: 1.2 m (4 ft)	1 cable	MMS: 866462; MP2: 865040; MP5: 865024, 865322; X2: 865039; MMX: 867036; Masimo® rainbow® Module: 867191; VS4: 863283; VS30: 863359, 863380; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304, 863322	•	•



989803186501



989803186531



Adapter cables

Masimo® LNCS® SpO₂ adapter cable

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803148221	n/a	Masimo® SET™ LNCS® adapter cable, adapts LNCS® sensors to Philips monitors Length: 3 m (10 ft)	1 cable	VS2: 863278, 863279; VS3: 863069, 863070, 863071, 863072, 863073, 863074; VS4: 863283; VSi: 863276, 863277; MMS: 862442; MP2: 865040; MP5: 865024, 865322, 865120; X2: 865039; X3: 867030; MX100: 867033; MMX: 867036; MX40: 865350, 865351; Pulse oximetry module: 862112; Masimo® SET™ Module: 867192; CM10: 863301, 863365; CM12: 863303, 863366; CM100: 863300; CM120: 863302; CM150: 863304, 863322; FM30: 862199; FM40: 865132; FM50: 865071	•	•



989803148221



Adapter cables

Masimo® RD SET® SpO₂ adapter cables

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206901	n/a	RD SET® MP-05 Philips Dual Keyed Cable; 5 ft.	1 cable	MX40: 865350, 865351, 865352; PM6000: 867311, 867313, 867315	•	•
989803206911	n/a	RD SET® MP-12 Philips Dual Keyed Cable; 12 ft.	1 cable		•	•

Adapter cables

Masimo® RD rainbow SET® SpO₂ adapter cables

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803206931	n/a	RD rainbow SET® 25R – 12 Cable; 12 ft.	1 cable	MX100: 867033; X3: 867030	•	•



989803206901



989803206911



989803206931



Continuous measurement
Reusable autoclavable probes

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803203561	n/a	Esophageal/rectal temperature probe (autoclavable), 12 FR, reusable, continuous monitoring, 2-prong plug Length: 3 m (10 ft)	1 probe	MP monitors: MP2: 865040; MP5: 865024; CM10: 863301, 863310; CM12: 863303, 863311; CM100: 863300, 863320; CM120: 863302, 863321; CM150: 863304, 863322; MMS: 862442; X2: 865039; X3: 867030; MX100: 867033; Temperature module: 862291; Hemodynamic ext.: M3012A; 867039 Capnography ext.: M3014A; 867040 Microstream™ ext.: 862393, 865377, 867041 MMX: 867036; Intrepid: 867172	•	•
989803203571	n/a	Esophageal/rectal temperature probe (autoclavable), 10 FR, reusable, continuous monitoring, 2-prong plug Length: 3 m (10 ft)	1 probe		•	•
989803203581	n/a	Skin surface temperature probe (autoclavable), tape-on, reusable, continuous monitoring, 2-prong plug Length: 3 m (10 ft), diameter: 9.5 mm (0.37 in)	1 probe		•	•



989803203561



989803203571



989803203581



Spot-check measurement

Predictive reusable temperature probes

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803143381	n/a	Reusable oral temperature probe and well kit, predictive temp monitoring 2.7 m (9 ft)	1 probe	VS2: 863278, 863081, 863082; VS3: 863070, 863073, 863074; VSi: 863276; VS30: 863359; 863380, MP5: 865024, 865322; VS4: 863283; PM6000: 867311, 867313, 867315	•	•
989803143391	n/a	Reusable rectal temperature probe and well kit, predictive temp monitoring 2.7 m (9 ft)	1 probe	VS2: 863278, 863081, 863082; VS3: 863070, 863073, 863074; MP5: 865024, 865322; VS4: 863283; PM6000: 867311, 867313, 867315	•	•



989803143381



989803143391



Patient monitoring batteries

Ordering code	Alternative ordering code	Image	Description	Shelflife	Packaging	Use with	CE	FDA
451261017451	M4607-60001	View >	Battery charger adapter		1 adaptor	X2: 865039; MP2: 865040	•	•
989803210521		View >	Rechargeable Lithium-ion battery, 10.7V 7.1Ah	4 months	1 battery	IntelliVue MP5: 865024, 865322, 865120; FM20: 862198; FM30: 862199; MX400: 866060; MX450: 866062; MX500: 866064; MX550: 866066	•	•
989803148701	M4607A	View >	Rechargeable lithium-ion battery, 10.8 volts	4 months	1 battery	X2: 865039; MP2: 865040	•	•
989803205811		View >	MX40 Lithium-ion Battery, 3 pack	2 months	3 batteries	MX40: 865350, 865351, 865352	•	•
989803174891		View >	MX40 battery adapter 3-pack, AA, disposable		1 adapter	MX40: 865350, 865351, 867146	•	•
989803205801		View >	MX40 Lithium-ion Battery, 1 pack	2 months	1 battery	MX40: 865350, 865351, 865352	•	•
989803174881		View >	VM lithium-ion battery	6 months	1 battery	VS <i>i</i> : 863275, 863276, 863277; VS2: 863278, 863279;	•	•
989803196521		View >	X3/MX100 lithium-ion battery	4 months	1 battery	X3: 867030; MX100: 867033	•	•
989803202201			X3/MX100 battery charger adapter		1 adapter		•	•
989803199221		View >	Lithium Ion Battery Pack 9 Cells	6 months	1 battery	860354, 860355, 860429, 860437, 863359, 863380, 863300, 863301, 863302, 863303, 863304, 863320, 863365, 863321, 863366, 863322, 860306, 860310, 860315, 860332, 860352, 860353, 863282, 863359, 863379, 863380	•	•
989803206871			Lithium Ion battery for PM6000 series, 6Ah 86.4Wh	5 months	1 battery	PM6000: 867311, 867313, 867315	•	



451261017451
M4607-60001



989803210521



989803148701
M4607A



989803205811



989803174891



989803205801



989803174881



989803196521



989803199221

Telemetry Pouches

Ordering code	Alternative ordering code	Description	Packaging	Use with	CE	FDA
989803210401		Protective Carry Pouch	50/Box	MX40, 865350, 865351, 865352, 865215	•	•
989803211031		Protective Carry Pouch	200/Box	MX40, 865350, 865351, 865352, 865215	•	•



989803210401
989803211031

Not all products are available in all geographies. Please check with your Philips representative to complete portfolio availability.



Patient monitors

Ordering code	Alternative ordering code	Description
862442	M3001A	IntelliVue multi-measurement server
863283		SureSigns VS4 vital signs monitor with NBP and SpO ₂
863300, 863320		Efficia CM100 patient monitor
863301		Efficia CM10 patient monitor
863365		Efficia CM10 patient monitor (Brazil)
863302,863321		Efficia CM120 patient monitor
863303		Efficia CM12 patient monitor
863366		Efficia CM12 patient monitor (Brazil)
863304, 863322		Efficia CM150 patient monitor
863359, 863380	M8105A	EarlyVue V530 vital signs monitor
865024		Intellivue MP5 portable patient monitor
865039		IntelliVue MMS X2 measurement module and patient monitor
865040		Intellivue MP2 portable patient monitor
865322		Intellivue MP5SC spot check patient monitor
865350		Intellivue MX40 wearable patient monitor (US only), 1.4 GHz Smart-hopping
865351, 867146		Intellivue MX40 wearable patient monitor, 2.4 GHz Smart-hopping
865352		Intellivue MX40 wearable patient monitor, 802.11a/b/g
866060		IntelliVue MX400 portable/bedside monitor
866062	M3001AL	Intellivue MX450 portable/bedside patient monitor
866064		Intellivue MX500 portable/bedside monitor
866066		Intellivue MX550 portable/bedside patient monitor
866462		IntelliVue multi-measurement server SLCP
867030		IntelliVue X3 measurement module and patient monitor
867033		Intellivue MX100 patient monitor



Patient monitors

Ordering code	Alternative ordering code	Description
867311		Patient Monitor 6100
867313		Patient Monitor 6300
867315		Patient Monitor 6500

Measurement modules

Ordering code	Alternative ordering code	Description
860335		XPER Flex Cardio medical equipment
860338		XPER Flex Cardio medical equipment
862119	M1026B	Anesthesia gas module
862187	M3014A	Capnography extension
862279	M1012A	Cardiac output module
866078	M1027B	IntelliVue EEG module
862277	M1006B	Invasive blood pressure module
862291	M1029A	Temperature module
862393	M3015A	Microstream™ CO, extension
862394	M3016A	Mainstream CO, 2nd press temp measurement extension
865021	M1014A	IntelliVue spirometry module
865377	M3015B	Microstream™ CO ₂ extension
865383		IntelliVue NMT module
866173		IntelliVue G7m anesthesia gas module
867039		Intellivue hemodynamic extension
867040		Capnography extension
867041		Microstream™ Extension
867191		Masimo rainbow' SET™ IV module
867409		FloTrac Module

Fetal monitors

Ordering code	Alternative ordering code	Description
862198	M2702A	Avalon FM20 fetal monitor
862199	M2703A	Avalon FM20 fetal monitor
865071	M2705A	Avalon FM50 fetal monitor
865132	M2704A	Avalon FM40 fetal monitor
867245		Avalon Toco+ MP transducer
867246		Avalon US transducer
867247		Avalon ECG/IUP transducer
867248		Avalon TOCO MP transducer
867249		Avalon Toco+ transducer
866075		Avalon Cableless Toco+ MP transducer
866076		Avalon Cableless ultrasound transducer
866077		Avalon Cableless ECG/IUP transducer



Defibrillators

Ordering code	Alternative ordering code	Description
866199		Efficia DFM100 defibrillator/monitor
867172		Heartstart Intrepid monitor/defibrillator

